

Analytical offers a unique new product, a 384-Well Desalting Plate for Peptides. This high-throughput peptide clean-up product works great for preparing peptide samples for mass spectrometry analysis.

Features and Advantages:

- 384-well plate with proprietary resin
- Works for both hydrophilic and hydrophobic peptides
- Removes urea and other salts from peptide samples with high recovery


High-throughput format allows simultaneous handling of 384 samples for peptide desalting with the use of centrifugation.



Cat. No.	Description	Qty
384PDP140	384-Well Peptide Desalting Plate	Each

Experimental Data:

HPLC chromatograms of peptide samples before (A) and after desalting (B) using 384-Well peptide desalting spin plate.



General Protocol:

STEP 1	Activate the hydrophobic, low carbon loading resin by adding 100μL of acetonitrile and centrifuge at 300 x g for 3 minutes and repeat the step.
STEP 2	Equilibrate the resin by adding 100μL of 0.1% trifluoroacetic acid (TFA) in water and centrifuge at 300 x g for 3 minutes and repeat the step for two additional times.
STEP 3	Add the sample (0.5μg – 30μg) in 50μL after diluting the sample 2-fold in 0.4% TFA in water and centrifuge at 300 x g for 3 minutes.
STEP 4	Wash the resin by adding 100 μL of 0.1% TFA in water and centrifuge at 300 x g for 3 minutes and repeat the step for two additional times.
STEP 5	Elute the bound peptides by adding 50 μL of 70% acetonitrile containing 0.1% TFA and centrifuge at 300 x g for 3 minutes and repeat the step.