

Flexi-Tier™

96-Well Collection Systems

- Unique design for maximum flexibility
- Adjusts to various heights
- Accommodates 0.5 - 2mL vials
- Compatible with all autosampler & robotic systems
- Can be centrifuged to 2500 RPMs
- Open bottom base allows immersion
- Total fluid contact with vials for even heat or cold transfer



Clear Flexi-Tier™

96-Position Flexi-Tier™ Complete Opaque Systems*

Heat Resistance up to 80°C

Cat. No.	Size	Material & Style	Base	Price
96FP05	0.5mL	Clear Conical Glass	Opaque	\$71
96FP10	1.0mL	Clear Conical Glass	Opaque	\$76
96FP15	1.5mL	Clear Conical Glass	Opaque	\$75
96F220	2.0mL	Clear Conical Glass	Opaque	\$87
96F220A	2.0mL	Amber Conical Glass	Opaque	\$87
96PV05	0.5mL	Conical Polypropylene	Opaque	\$51
96PV10	1.0mL	Conical Polypropylene	Opaque	\$57
96PV15	1.5mL	Conical Polypropylene	Opaque	\$67
96PV20	2.0mL	Conical Polypropylene	Opaque	\$64
96BB05	0.5mL	Clear Flat Bottom Glass	Opaque	\$71
96BB10	1.0mL	Clear Flat Bottom Glass	Opaque	\$76
96BB15	1.5mL	Clear Flat Bottom Glass	Opaque	\$64
96FP20	2.0mL	Clear Flat Bottom Glass	Opaque	\$82



Opaque Flexi-Tier™

96-Position Flexi-Tier™ Complete Clear Systems*


Heat Resistance up to 150°C

Cat. No.	Size	Material & Style	Base	Price
96FP05-C	0.5mL	Clear Conical Glass	Clear	\$73
96FP10-C	1.0mL	Clear Conical Glass	Clear	\$76
96FP15-C	1.5mL	Clear Conical Glass	Clear	\$75
96F220-C	2.0mL	Clear Conical Glass	Clear	\$98
96F220A-C	2.0mL	Amber Conical Glass	Clear	\$88
96PV05-C	0.5mL	Conical Polypropylene	Clear	\$51
96PV10-C	1.0mL	Conical Polypropylene	Clear	\$57
96PV15-C	1.5mL	Conical Polypropylene	Clear	\$67
96PV20-C	2.0mL	Conical Polypropylene	Clear	\$64
96BB05-C	0.5mL	Clear Flat Bottom Glass	Clear	\$65
96BB10-C	1.0mL	Clear Flat Bottom Glass	Clear	\$64
96BB15-C	1.5mL	Clear Flat Bottom Glass	Clear	\$64
96FP20-C	2.0mL	Clear Flat Bottom Glass	Clear	\$76


Compatible with all autosampler
& robotic systems



* All Systems Include Pre-Scored PTFE/Silicone Flat Liner (96FSC1) and Cover.




Conical Bottom Clear Glass Vials



Conical Bottom Amber Glass Vials



Flat Bottom Clear Glass Vials



Flat Bottom Amber Glass Vials

Glass Inserts in Loader

Cat. No.	Type	Style	Volume	Qty	Price
96FV05	Conical	Clear	0.5mL	96	\$41
96FV10	Conical	Clear	1.0mL	96	\$42
96FV15	Conical	Clear	1.5mL	96	\$41
96FV20	Conical	Clear	2.0mL	96	\$47
96FV05A	Conical	Amber	0.5mL	96	\$44
96FV10A	Conical	Amber	1.0mL	96	\$44
96FV15A	Conical	Amber	1.5mL	96	\$57
96FV20A	Conical	Amber	2.0mL	96	\$47
96VL05	Flat Bottom	Clear	0.5mL	96	\$20
96VL10	Flat Bottom	Clear	1.0mL	96	\$20
96VL15	Flat Bottom	Clear	1.5mL	96	\$20
96VL20	Flat Bottom	Clear	2.0mL	96	\$22
96VL05A	Flat Bottom	Amber	0.5mL	96	\$38
96VL10A	Flat Bottom	Amber	1.0mL	96	\$28
96VL15A	Flat Bottom	Amber	1.5mL	96	\$28
96VL20A	Flat Bottom	Amber	2.0mL	96	\$28

Silanized Glass Inserts in Loader

96FV05-S	Conical	Clear	0.5mL	96	\$54
96FV10-S	Conical	Clear	1.0mL	96	\$55
96FV15-S	Conical	Clear	1.5mL	96	\$54
96FV20-S	Conical	Clear	2.0mL	96	\$59
96FV20A-S	Conical	Amber	2.0mL	96	\$60
96VL05-S	Flat Bottom	Clear	0.5mL	96	\$35
96VL10-S	Flat Bottom	Clear	1.0mL	96	\$35
96VL15-S	Flat Bottom	Clear	1.5mL	96	\$35
96VL20-S	Flat Bottom	Clear	2.0mL	96	\$37
96VL20A-S	Flat Bottom	Amber	2.0mL	96	\$54

Bulk, Loose Glass Inserts

Cat. No.	Type	Style	Volume	Qty	Price
96F005	Conical	Clear	0.5mL	100	\$34
96F010	Conical	Clear	1.0mL	100	\$32
96F015	Conical	Clear	1.5mL	100	\$34
96F020	Conical	Clear	2.0mL	100	\$39
96F005A	Conical	Amber	0.5mL	100	\$37
96F010A	Conical	Amber	1.0mL	100	\$37
96F015A	Conical	Amber	1.5mL	100	\$37
96F020A	Conical	Amber	2.0mL	100	\$40
96F605	Flat Bottom	Clear	0.5mL	100	\$13
96F610	Flat Bottom	Clear	1.0mL	100	\$13
96F615	Flat Bottom	Clear	1.5mL	100	\$14
96F620	Flat Bottom	Clear	2.0mL	100	\$15
96F625	Flat Bottom	Clear	2.5mL	100	\$18
96F605A	Flat Bottom	Amber	0.5mL	100	\$18
96F610A	Flat Bottom	Amber	1.0mL	100	\$19
96F615A	Flat Bottom	Amber	1.5mL	100	\$17
96F620A	Flat Bottom	Amber	2.0mL	100	\$21

Polypropylene Inserts in Loader

Cat. No.	Type	Style	Volume	Qty	Price
96PL05	Conical	PP	0.5mL	96	\$22
96PL10	Conical	PP	1.0mL	96	\$22
96PL15	Conical	PP	1.5mL	96	\$24
96PL20	Conical	PP	2.0mL	96	\$31

PP Inserts w/ Siliguard in Loader

96PL05GC	Conical	PP	0.5mL	96	\$35
96PL10GC	Conical	PP	1.0mL	96	\$36
96PL15GC	Conical	PP	1.5mL	96	\$37
96PL20GC	Conical	PP	2.0mL	96	\$38

Bulk, Loose PP Inserts

Cat. No.	Type	Style	Volume	Qty	Price
96P005	Conical	PP	0.5mL	100	\$15
96P010	Conical	PP	1.0mL	100	\$16
96P015	Conical	PP	1.5mL	100	\$18
96P020	Conical	PP	2.0mL	100	\$20

Bulk, Loose PP Inserts w/ Siliguard


96P005GC	Conical	PP	0.5mL	100	\$29
96P010GC	Conical	PP	1.0mL	100	\$30
96P015GC	Conical	PP	1.5mL	100	\$31
96P020GC	Conical	PP	2.0mL	100	\$34




Conical Bottom Polypropylene Vials

Flexi-Tier™ Collection Plate Accessories

Cat. No.	Description	Qty	Price
96FC00	Opaque Cover	5	\$34
96FS00	Opaque Spacer	1	\$10
96FB00	Opaque Base Plate	1	\$13
96RB00	Opaque Round Base Plate	1	\$55
96FS00-C	Clear Spacer	1	\$10
96FC00-C	Clear Cover	5	\$34
96FB00-C	Clear Base Plate	1	\$13



96FB00 Opaque Open Bottom Base Plate





96RB00 Opaque Open Round Bottom Base Plate



Plug Caps for 1mL Polypropylene Tubes

- For mixing and long-term storage
- Not recommended for autosampler

Cat. No.	Description	Qty	Price
96812-2	Strips of 12 Plug Caps	80	\$17
96128-2	Strips of 8 Plug Caps	120	\$13
96962-2	Individual Plug Caps	960	\$50




Cap Mats For Flexi-Tier & Aluminum Collection Systems

Cat. No.	Description	Qty	Price
965007	Ultra Thin Round Well, Soft Silicone/PTFE, Purple	5	\$35
963675	Round Well, Silicone/PTFE, Green	5	\$65
96057	1mL Autosampler Compatible Round Cap Mat, Polypropylene with 'X' on Surface, Clear	20	\$54
967702	Round Well Molded Soft Silicone, Blue	5	\$62
964007	Ultra Thin Round Well w/ Molded Soft Silicone, Clear	5	\$35
96FSC2	Molded PTFE/Silicone Liner, White	5	\$148
96FSC3	Pre-Scored Molded PTFE/Silicone Liner, White	5	\$133
96FSC4	Molded Silicone Liner (Standard with Aluminum Systems)	5	\$84
96FSC5	Pre-Scored Molded Silicone Liner, White	5	\$87
96FSC3STRIP	Pre-Scored Molded PTFE/Silicone 8 Position Strip Liner	12	\$47
96FSC3PLUG	Pre-Scored Molded PTFE/Silicone Single Plug	100	\$55
96FSC2STRIP	Molded PTFE/Silicone 8 Position Strip Liner, White	12	\$45
96FSC2-12	Molded PTFE/Silicone 12 Position Strip Liner, White	8	\$44
96FSC2PLUG	Molded PTFE/Silicone Single Plug, White	100	\$53



Advantage™ 96-Well Standard Aluminum Base Plates

Cat. No.	Description	Price
96AB11	Aluminum Base Plate Slotted/Open Bottom	\$111
96AB22	Aluminum Base Plate Solid Bottom	\$115
96AB74	Aluminum Open Round Bottom Base	\$113



96AB11 Aluminum Slotted Base Plate

Advantage™ 96-Well Tall Aluminum Base Plate

- Superb conductivity for uniform vial heating and cooling, and can be immersed in heating/ cooling liquids for total fluid contact with the vials.
- System fits all liquid handling systems, including most autosamplers.

Cat. No.	Description	Price
96A002	Tall Aluminum Base Plate	\$854



Tall Aluminum Base Plate – provides more consistent heating/cooling of 1.5mL and 2mL vials.

Snap/Crimp Vials and Caps for Flexi-Tier™ Systems



Snap & Seal Vials for Flexi-Tier System

Cat. no.	Description	Qty	Price
96922	Conical Clear Glass 0.5mL Snap & Seal Vial	100	\$45
96936	Conical Clear Glass 0.75mL Snap & Seal Vial	100	\$45
96938	Tapered Clear Glass 0.75mL Snap & Seal Vial	100	\$49
96946	Conical Clear Glass 1.0mL Snap & Seal Vial	100	\$45
96948	Tapered Clear Glass 1.0mL Snap & Seal Vial	100	\$49

Snap & Seal Caps for Flexi-Tier Conical and Tapered Vials

Cat. No.	Description	Qty	Price
96030	PTFE Lined Cap	100	\$13
96040	PTFE/Butyl Rubber Lined Cap	100	\$16
96080	PTFE/Silicone Lined Cap	100	\$22
96060	PTFE/Silicone/PTFE Lined Cap	100	\$28
96070	Pre-Slit PTFE/Silicone Lined Cap	100	\$29

Snap & Seal Vials for Flexi-Tier System – Silanized

Cat. no.	Description	Qty	Price
96922S	Silanized Conical Clear Glass 0.5mL Snap & Seal Vial	100	\$59
96936S	Silanized Conical Clear Glass 1.0mL Snap & Seal Vial	100	\$59
96938S	Silanized Tapered Clear Glass 1.0mL Snap & Seal Vial	100	\$63
96946S	Silanized Conical Clear Glass 1.5mL Snap & Seal Vial	100	\$59
96948S	Silanized Tapered Clear Glass 1.5mL Snap & Seal Vial	100	\$63

Flexi-Tier™ Aluminum Collection Systems

Advantage™ Aluminum 96-Well Collection System for High Throughput

- Solid aluminum base design is perfect for platform heating/cooling devices.
- Superb conductivity for uniform vial heating and cooling, and can be immersed in heating/cooling liquids for total fluid contact with the vials.
- Ideal for volatile solvents, salt screening, high throughput solubility & physical characterization.
- Samples are securely contained without possibility of crosstalk.
- System fits all liquid handling systems, including most autosamplers.

Cat. No.	Description	Qty	Price
96A105	Complete 0.5mL Aluminum Block System with Conical Clear Glass Inserts	Each	\$248
96A110	Complete 1.0mL Aluminum Block System with Conical Clear Glass Inserts	Each	\$249
96A115	Complete 1.5mL Aluminum Block System with Conical Clear Glass Inserts	Each	\$259
96A120	Complete 2.0mL Aluminum Block System with Conical Clear Glass Inserts	Each	\$265
96F105	Complete 0.5mL Aluminum Block System with Flat Bottom Clear Glass Inserts	Each	\$231
96F110	Complete 1.0mL Aluminum Block System with Flat Bottom Clear Glass Inserts	Each	\$231
96F115	Complete 1.5mL Aluminum Block System with Flat Bottom Clear Glass Inserts	Each	\$231
96F120	Complete 2.0mL Aluminum Block System with Flat Bottom Clear Glass Inserts	Each	\$231

Replacement Parts

96A001	Aluminum 96-Well Solid Bottom Base Plate w/ Tapped Holes in Each Corner	Each	\$120
96AB77	Aluminum 96-Well Round Hole Bottom Base Plate, Tapped Holes in Each Corner	Each	\$115
96B002	Aluminum 96-Well Block System Slotted Cover	Each	\$70
96050	Stainless Steel Screws for 0.5mL Aluminum 96-Well Block System	4	\$6
96010	Stainless Steel Screws for 1.0mL Aluminum 96-Well Block System	4	\$6
96015	Stainless Steel Screws for 1.5mL Aluminum 96-Well Block System	4	\$6
96020	Stainless Steel Screws for 2.0mL Aluminum 96-Well Block System	4	\$6

