

384-Well Collection Plates

Our Advantage™ 384-Well 1536-Well collection plates and accessories are designed for high throughput analysis and are available in an assortment of volume sizes and colors. The ability to color code allows you to quickly and easily identify the various stages of your samples whether they have been processed through protein precipitation, de-salting or dilutions.

Our polypropylene plates are chemical and solvent resistant and are fully compatible with automated and robotic systems.


384-Well Polypropylene Collection Plates

Cat. No.	Description	Qty/Pkg
38435	0.35mL, Clear PP, 58µL per well, Square Well, Spherical Bottom	60
38120	120µL, Clear Polypropylene, Square Well, Flat Bottom	24
38200	200µL, Polypropylene, Clear, Square Well, Flat Bottom	24
38200B	200µL, Polypropylene, Blue, Square Well, Flat Bottom	24
38200G	200µL, Polypropylene, Green, Square Well, Flat Bottom	24
38200R	200µL, Polypropylene, Red, Square Well, Flat Bottom	24
38200Y	200µL, Polypropylene, Yellow, Square Well, Flat Bottom	24
38400	400µL, Polypropylene, Clear, Square Well, Flat Bottom	24
38400B	400µL, Polypropylene, Blue, Square Well, Flat Bottom	24
	400µL, Polypropylene, Green, Square Well, Flat Bottom	24
38400R	400µL, Polypropylene, Red, Square Well, Flat Bottom	24
38400Y	400µL, Polypropylene, Yellow, Square Well, Flat Bottom	24

384280	12-130µL Working Volume, PP, Clear, Square Well, V-Bottom	100
384241	1-90µL Working Volume, PP, Clear, Square Well, V-Bottom	100
384201	15-145µL Working Volume, PP, Clear, Square Well, Flat Bottom	100

384 Square Well Purple Pattern Adhesive PTFE (Teflon®) Sealing Film


- Ultra Thin, 2mil
- Superior Solvent Compatibility
- Higher Chemical Resistance
- Greater Tackiness
- DMSO Resistant




Cat. No.	Description	Qty
38402	384 Square Well Pattern Adhesive PTFE Sealing Film, Purple	100
38402SL	384 Square Well Pattern Adhesive PTFE Sealing Film, Pre-Slit, Red	100

1536-Well Collection Plates – Square Well

Cat. No.	Description	Qty
15020	Polystyrene, Square Top/Square Bottom, 20µL	60





Advantage™ 384-Well Collection Plates – Round Well

Cat. No.	Description	Qty
38650	65µL, Polystyrene, Clear , Round Well, Flat Bottom	100
38650B	65µL, Polystyrene, Black , Round Well, Flat Bottom	100
38650W	65µL, Polystyrene, White , Round Well, Flat Bottom	100
38450	65µL, Polystyrene, Clear , Round Well, V-Bottom	100
38450B	65µL, Polystyrene, Black , Round Well, V-Bottom	100
38450W	65µL, Polystyrene, White , Round Well, V-Bottom	100
38465C	65µL, Polypropylene, Clear , Round Well, Flat Bottom	100
38465B	65µL, Polypropylene, Black , Round Well, Flat Bottom	100
38465W	65µL, Polypropylene, White , Round Well, Flat Bottom	100

Cap Mats & Accessories

384-Well Cap Mats

Cat. No.	Description	Qty
96384	Molded 384-Well Soft Silicone Square Well Cap Mat	5
96384SL	Molded 384-Well Pre-Slit Soft Silicone Square Well Cap Mat	5

Standard Sealing Films


- Temperature range -40°C to 120°C
- Sterile films also available - call for pricing

Cat. No.	Description	Qty
968120	Aluminum Foil Film for PCR, 1.5mils, Pierceable	100
964012	Clear Polypropylene Film for PCR, 2mils, <i>not pierceable</i>	100


Rubber Roller

- Secures cap mat to plate



Cat. No.	Description	Qty
96400	4" Rubber Roller	1



96384SL



968120 (foil)



Light Up Your 384-well Collection Plate Samples

LightOne Pro Illuminator

- Quickly identify and select the area you want to pipette
- Customizable: highlight a single well or an entire row or column
- Foot pedal control for easy, hands-free navigation
- Use USB port for downloading templates and customized files for cherry picking
- 384-Well or 96-Well

Cat. No.	Description
LOI-384	LightOne Pro

For pricing, see analytical-sales.com