



analytical



48 Vial Cap Mat



2026



High Throughput Screening for High Throughput Sampling

NEW! 48-position Silicone/PTFE Cap Mats

For 9mm, 12 x 32mm HPLC Vials in the **Waters® Acquity® Autosampler Tray**

- *Efficiently Cap 48 HPLC Vials at Once!*
- *Great for Dissolution Assays*

Made in the USA

Cat. No.	Description	Qty
99970	Blue , Silicone/PTFE Cap Mats for 2mL, 12x32mm HPLC Vials	5
99975	Clear, Pre-Slit , Silicone/PTFE Cap Mats for 2mL, 12x32mm HPLC Vials	5

For pricing, see analytical-sales.com