



20mm Crimp Style Head Space Vials

Crimp style vials, made of borosilicate glass

For continuous results it is very important the head space vial and cap form a tight seal. In this type of analysis, the sample is sealed and heated to very high temperatures so the concentration in the gas phase (head space) is equal to the amount in the liquid phase. An aliquot from the headspace is then analyzed.



62700

Cat. No.	Description	Qty
62321-3	6mL Clear Glass Crimp Vial, Bevel Top / Flat Bottom, 22 x 38mm	100
62101-3	10mL Clear Glass Crimp Vial, Bevel Top / Flat Bottom, 23 x 46mm	100
62201-3	20mL Clear Glass Crimp Vial, Bevel Top / Flat Bottom, 23 x 75mm	100

Recommended for PE Autosamplers

67220-3	20mL Amber Glass Crimp Vial, Bevel Top / Round Bottom, 23 x 75mm	100
---------	---	-----

Recommended for CTC Autosamplers

62282-3	6mL Clear Glass Crimp Vial, Beveled Top / Round Bottom, 22 x 38mm	100
61182-3	10mL Amber Glass Crimp Vial, Beveled Top / Round Bottom, 23 x 46mm	100
61182-3-Case	10mL Amber Glass Crimp Vial, Beveled Top / Round Bottom, 23 x 46mm	1000



62800

20mm Crimp Caps with Liners

Available in Aluminum, Aluminum/Tin and Magnetic Caps

Magnetic Caps are made from a special polished steel alloy for optimum reliability.

Cat. No.	Description	Qty
62700	20mm Aluminum Cap, Ultra Low Bleed, Beige PTFE / White Silicone	100
62800	20mm Blue Magnetic Aluminum / Tin Cap, Ultra Low Bleed, Beige PTFE / White Silicone	100
62900	20mm Magnetic Aluminum Cap, Ultra Low Bleed, Beige PTFE / White Silicone	100

For pricing, see analytical-sales.com

See Screw Thread Head Space Vials on reverse




18mm Screw Thread Head Space Vials

Cat. No.	Description	Qty
62546	10mL Clear Glass Flat Bottom Vial, 22.5 x 46.0mm	100
62546A	10mL Amber Glass Flat Bottom Vial, 22.5 x 46.0mm	100
62546S	10mL Clear Glass Flat Bottom Vial, 22.5 x 46.0mm - Silanized	100
62546S-CASE	10mL Clear Glass Flat Bottom Vial, 22.5 x 46.0mm - Silanized	1000
62346A	10mL Amber Glass Round Bottom Vial, 22.5 x 46.0mm	100
72546	20mL Clear Glass Flat Bottom Vial, 22.5 x 75.5mm	100
72546A	20mL Amber Glass Flat Bottom Vial, 22.5 x 75.5mm	100
72546S	20mL Clear Glass Flat Bottom Vial, 22.5 x 75.5mm - Silanized	100
72546S-CASE	20mL Clear Glass Flat Bottom Vial, 22.5 x 75.5mm - Silanized	1000




18mm Magnetic Screw Thread Closures with Ultra Low Bleed Silicone Liners

Cat. No.	Description	Qty
62515	18mm Magnetic Screw Cap with Teflon Liner, 0.010"	100
62518	18mm Magnetic Screw Cap with Blue Teflon / White Silicone, 0.080"	100
62517	18mm Magnetic Screw Cap with Blue Teflon / White Silicone, 0.080"	1000



Custom Head Space Vial Holders

Cat. No.	Description	Qty
6233	Head Space Vial Holders for 18mm & 20mm Vials, 3 x 3	1
6244	Head Space Vial Holders for 18mm & 20mm Vials, 4 x 4	1
6255	Head Space Vial Holders for 18mm & 20mm Vials, 5 x 5	1
6266	Head Space Vial Holders for 18mm & 20mm Vials, 6 x 6	1



Crimpers & Decappers

Cat. No.	Description	Qty
90200	Vial Crimper for 20mm OD Vials	1
90100	Decapping Pliers for All Vials	1

For pricing, see analytical-sales.com