



627350



627400



627550



627650

Stainless Steel Vial Racks

24, 48 and 96 well formats

SS Vial Racks enable workflows or situations that require vials to be exposed to open air or submerged in liquid. For example:

- Sterilization of vials and holder by autoclave
- Sonicator baths
- Incubators
- Long term storage (refrigerator/freezer)

Features

- High-quality Stainless Steel
- Automation compatible (SLAS footprint)
- Open sides provide maximized exposure for gases or liquids
- Can be cleaned and reused, making it a “green” alternative to disposable plastic trays/holders
- Vials can be sealed with silicone cap mats, TPE plugs, appropriate screw/crimp caps, etc...

Cat. No.	Description	Qty
627350	96 Position Stainless Steel Vial Rack for 8mm OD Vials	Each
627400	48 Position Stainless Steel Vial Rack for 12mm OD Vials - VHP Compatible	Each
627550	24 Position Stainless Steel Vial Rack for 1 Dram (15mm OD) Vials (18mm well spacing)	Each
627650	24 Position Stainless Steel Vial Rack for 2 Dram (17mm OD) Vials (18mm well spacing)	Each