

Storage Boxes for 3-5mL Storage Tubes

- Made of high impact polystyrene
- Temperature range: -80°C to 80°C. Not autoclaveable.

Cat. No.	Description	Qty
96581B	Storage Box, Holds 81 (3-5mL) Sample Tubes, Blue	10



100 Well Microtube Storage Boxes For Microtubes & Cryogenic Vials

- Autoclavable
- High Transparency

Cat. No.	Description	Qty
90006	Plastic Storage Boxes Hold 1.5/2.0 mL Cryovials or Microtubes	5

Silicone Lab Mats

The Lab Mat is an economical solution to help keep benchtops clean and safe from stains, spills and wear

- Large working surface area accomodates benchtop equipment, weigh boats, stir bars, scales, etc.
- Protect countertops from hot items up to 200°C (392°F)
- Reversible, with two color/design options available


Silicone Lab Mats

Cat. No.	Description	Qty
LABMATP	Silicone Lab Mat, Purple/Gray	Each
LABMATY	Silicone Lab Mat, Yellow/Gray/Blue	Each



• Lip design around all edges enables liquid retention of approximately 300ml of fluid

• Internal metal inserts for retention of magnetic objects, such as stir bars to help keep items organized



Combi-Rac™ Versatile Sample Rack Systems

Uses minimal counter space and offers maximum flexibility. Made of autoclavable polypropylene. One side holds eighty 10 & 12mmOD test tubes or HPLC vials. Flip the tray over and it holds sixty micro-tubes from 0.5 to 2mL volumes. Overall size measures 9"L x 2.5"W x 2.5"H. Comes with 2 removable handles and is available in a variety of colors.


Cat. No.	Description	Color	Qty
8B165	Combi-Rac™	Blue	10
8G165	Combi-Rac™	Green	10
8R165	Combi-Rac™	Red	10
8P165	Combi-Rac™	Pink	10
8O165	Combi-Rac™	Orange	10
8Y165	Combi-Rac™	Yellow	10
8C165	Combi-Rac™	Assorted	10
8L165	Combi-Rac™ Lid	Transparent	10



Mini Combi-Rac™

- 3" x 3" autoclavable PP rack
- Flip it – one side holds 25 0.5-2.0mL tubes or vials, other side holds 16 0.2-0.5mL micro centrifuge tubes

Cat. No.	Description	Color	Qty
5B025	Mini Combi-Rac™	Blue	10
5R025	Mini Combi-Rac™	Red	10
5Y025	Mini Combi-Rac™	Yellow	10



4-Way Rack

- Link together to cater for your application needs
- Create as large a rack as you require for your bench top applications
- Made of polypropylene
- Autoclavable
- Assorted colors
- 6.85" x 3.74" x 2.04" (17.4 x 9.5 x 5.2 cm)

Cat. No.	Description	Tubes	Wells	Well Dia.	Qty
4WRack	4-Way Rack	50mL	4	30mm	5
		15mL	12	17.5mm	
		0.5mL	32	8.3mm	
		1.5/2.0mL	32	11.6mm	



Tube Racks

Cat. No.	Description	Color	Qty
813600B	Racks for Tubes up to 13mm	Blue	10
813600G	Racks for Tubes up to 13mm	Green	10
816300B	Racks for Tubes up to 16mm	Blue	10
816300G	Racks for Tubes up to 16mm	Green	10
825300B	Racks for Tubes up to 25-30mm	Blue	10
825300G	Racks for Tubes up to 25-30mm	Green	10

One-Well 50 mL Connecting Tube Racks

- Interlocking - form a line of 50mL tubes, or use as a single, stand-alone rack
- Molded circular base ring keeps conical or round-bottom tubes in place
- Autoclavable



Cat. No.	Description	L x W x H (mm)	Colors	Qty
CTR1050	1-Well 50 mL Tube Racks	83 x 70 x 67 mm	Assorted	5



One-Well 50 mL Single Tube Racks

- Holds conical- or round-bottom tubes.
- Autoclavable, epoxy-coated steel
- Base ring keeps tubes in place
- Measures 90 x 67 x 67 mm

Cat. No.	Description	Color	Qty
WR1050	Epoxy-Coated Wire Rack	Yellow	1

96-Well PCR® Racks with Hinged Lid

Polypropylene racks hold 96 (0.2mL) single tubes or 8/12-tube strips. Great for pre or post PCR® sampling applications. Imprinted numbers and letters for easy identification. Removable hinged lid, available in assorted colors.

PCR® is a registered trademark of Hoffman-LaRoche Inc.




Cat. No.	Description	Color	Qty
PCRGR	96-Well PCR® Racks	Green	5
PCRBL	96-Well PCR® Racks	Blue	5
PCRYL	96-Well PCR® Racks	Yellow	5
PCRNA	96-Well PCR® Racks	Natural	5
PCRAS	96-Well PCR® Racks	Assorted	5



Flatpack Freezer Boxes

Reusable, economical and moisture resistant plastic

- Made of polypropylene material
- Easy self-assemble boxes shipped flat
- More durable than cardboard boxes as they are not susceptible to moisture
- Significantly reduces the chance of mold
- Dividers included with the boxes



Freezer boxes are supplied "flat-packed" eliminating the hassle of stacking and reduces the amount of space needed for storage.

Thin Polypropylene Film 1.5-2 mL

• Can be used in mechanical freezers and vessels down to 196°C (vapor phase only)

Cat. No.	Tube	Wells	Well Dia.	L x W x H	Color	Qty
FFB81W	1.5/2.0 mL	81	12.1 mm	4.95" x 4.92" x 1.92" (12.5 x 12.4 x 4.9 cm)	White	10
FFB81B					Blue	
FFB81G					Green	
FFB81P					Purple	



Corrugated Polypropylene Box with Dividers

• Store down to temperatures of -80°C

Cat. No.	Tube	Wells	Well Dia.	LxWxH (mm)	Color	Qty
FFB81	0.5 mL	81	12 mm	133 x 133 x 51	Natural	10

Cardboard Micro Tube Boxes

Safe micro tube storage

These low profile boxes are available in two sizes to accommodate standard 0.5mL and 0.2mL micro tubes and vessels of a similar size. They come complete with well dividers in packs of 10.



- Store in mechanical freezers or vapor phase LN2*
- Will fit into most standard freezer racking
- 0.2 mL tube (9.5 mm max) version has 169 wells in a 13 x 13 format
- 0.5 mL tube (12 mm max) version has 100 wells in a 10 x 10 format

Cat. No.	Description	Tubes	Wells	LxWxH (mm)	Qty
MTB169	Cardboard Micro Tube Boxes	0.2 mL	169	133 x 133 x 29	10
MTB100	Cardboard MicroTube Boxes	0.5 mL	100	133 x 133 x 41	10

Storage Racks for 80-Well Racks


Strong, clear acrylic construction

For use at room temperature only - do not freeze.

- Holds up to eight 80-well microtube racks filled with either 1.5 or 2.0 mL microtubes



Cat. No.	Description	L x W x H (mm)	Qty
STRV80	Vertical Rack	157 x 236 x 234	1
STRH80	Horizontal Rack	310 x 236 x 122	1



Benchtop Dispensing Bins

Large, angled openings with oversized lip on lid for easy access. Green bin includes a removable divider. Non-skid rubber feet and pre-drilled holes for wall mounting (hardware included).

Cat. No.	Description	Color	Qty
28322	180W x 155H x 165mmD, Dual Compartment	Neon Green	1
28321	155W x 155H x 170mmD, Single Compartment	Neon Blue	1

Crystal-clear! See your small item inventory!



Workstation Dispensing Bins

Flip-top covers make filling and dispensing easy and keep products dust-free. With non-skid rubber feet.

Cat. No.	Description	Qty
97321	184W x 194H x 232mm D, One Bin	1
97322	150W x 177H x 368mm D, Two Bins	1
97323	217W x 170H x 217mm D, Three Bins	1
97324	287W x 170H x 217mm D, Four Bins	1

Pour-boat Weighing Dishes

Pour-boat Weighing Dishes

- Pour-spout design easily dispenses samples without flexing the dish.
- Flat-bottom dishes are made from 12-mil thick polystyrene for better strength, stability, and handling



Cat. No.	Description	Size	L x W x H	Qty
PBWS	Pour-boat Weighing Dish	Small	43 x 58 x 13 mm	500
PBWM	Pour-boat Weighing Dish	Medium	83 x 132 x 26 mm	500
PBWL	Pour-boat Weighing Dish	Large	108 x 183 x 26 mm	500

Diamond Weighing Boats

- Diamond shape that easily flexes for accurate pouring
- Wide, flat bottom to resist tipping



Cat. No.	Description	Size	L x W x H	Qty
DWB5	Diamond Weighing Boats	5 mL	55 x 35 x 6 mm	500
DWB25	Diamond Weighing Boats	25 mL	71 x 46 x 14 mm	500
DWB30	Diamond Weighing Boats	30 mL	78 x 56 x 14 mm	500
DWB100	Diamond Weighing Boats	100 mL	119 x 90 x 19 mm	500

Parafilm® Sealing Film

Colorless, moisture-proof film molds to surface

Parafilm® M is a semi-transparent, flexible, thermoplastic film with two-way stretch.

- Non-toxic and free from plasticizers, primarily composed of polyolefins and paraffin waxes offering good resistance against acids, alkalis and some organic solvents – not recommended for use with acetone, diethyl ether or tetrachloride
- Will stretch 200% and adhere to irregular shapes and surfaces
- Disposable



Cat. No.	Description	Film Width	Roll Length
PRS5075	Parafilm® Sealing Film	50 mm (2-inch)	75m
PRS10038	Parafilm® Sealing Film	100 mm (4-inch)	38m
PRS10075	Parafilm® Sealing Film	100 mm (4-inch)	75m

Parafilm® Dispenser

Heavy-duty acrylic for increased durability

Store, dispense, and cut Parafilm®, tape, labels, and Tough-Tags®.

- Strong, crystal-clear acrylic dispenser makes your film supply apparent and features a serrated cutting blade
- Dispenser holds one or two 50 mm (2-inch) rolls, or one 100 mm (4-inch) roll




Cat. No.	Description	w x h x d (mm)
PRSDISP	Acrylic Dispenser	119 x 168 x 175

Basins & Trough Plates

Multi-Channel Disposable Solution Basins

Constructed of white polystyrene or clear PVC, these disposable basins have angled sides and V-shaped bottoms allowing nearly complete sample withdrawal. 55mL volume.

Cat. No.	Description	Qty
DSBPVC	Clear, PVC Disposable Basins	100
DSBWPS	White, Polystyrene Disposable Basins	100



Advantage 8 and 12 Channel V-Troughs

Uniquely designed to minimize dead volume while you add standards, wash, rinse or use reagents in your high throughput clean ups. Can be used in liquid handling systems like TomTec etc.

Cat. No.	Description	Qty
96012	8 Row V-Trough, White Polystyrene	5
96013	12 Column V-Trough, Clear PETG	5
96FB00	Opaque Base Plate for 96012 Tray	Each

For added support and stability, try using a base plate (96FB00) under your V-troughs (96012 and 96013).

