




BLACK TrueTaper™

For Light Sensitive Compounds



TrueTaper™ 96-Well Collection Plates were designed for chromatographers faced with limited sample quantity. TrueTaper™ incorporates a smooth, tapered well shape which allows the entire sample to gravitate to the bottom of the well.

- Square wells with unique **TAPERED** round bottoms for maximum sample recovery
- 100µL tapered bottom volume
- 2.0mL total volume capacity

Cat. No.	Description	Qty
968820BL	TrueTaper 96-Well Collection Plate, Black	25

- Also available with **Siliguard** non-binding coating -

For pricing, see analytical-sales.com