

**Canary Safe™****COLOR-Coded  
Mobile Phase Safety Caps*****Easily identify your solvent lines!***

Here's just one example of our color-coded caps assigned to solvent lines:

- Red cap = Methanol
- Mustard cap = Acid
- Blue cap = Water
- Green cap = Base
- Purple cap = Acetonitrile

Of course, you can come up with your own color system to fit your needs!

**Color-Coded Mobile Phase  
Bottle Safety Caps****GL 45 Thread**

Cat. No.	Description	Color
CC1001	1 Standard Port for 1/8" OD Tubing	Purple
CC1002	2 Standard Ports for 1/8" OD Tubing	Purple
CC1001R	1 Standard Port for 1/8" OD Tubing	Red
CC1002R	2 Standard Ports for 1/8" OD Tubing	Red
CC1001M	1 Standard Port for 1/8" OD Tubing	Mustard
CC1002M	2 Standard Ports for 1/8" OD Tubing	Mustard
CC1001B	1 Standard Port for 1/8" OD Tubing	Dark Blue
CC1002B	2 Standard Ports for 1/8" OD Tubing	Dark Blue
CC1001G	1 Standard Port for 1/8" OD Tubing	Green
CC1002G	2 Standard Ports for 1/8" OD Tubing	Green
CX0001	Replacement Filter Check Valve for Mobile Phase Safety Caps	

**3 and 4 port caps are also available**

CX0001