

GL38 Caps



GL38 Mobile Phase Bottle Safety Caps



- Screws directly onto GL38 threaded amber solvent bottles (*NOTE: Most 4L solvent bottles are GL38 threads*)
- No thread adapter needed!
- Comes with Air Check Valve and Fittings

GL38 Mobile Phase Bottle Safety Caps

Cat. No.	Description
CCGL3801	Canary-Safe Cap, GL38, 1 Standard Tubing Port
CCGL3802	Canary-Safe Cap, GL38, 2 Standard Tubing Ports
CCGL3803	Canary-Safe Cap, GL38, 3 Standard Tubing Ports



Replacement Filter Check Valve

- Fits Canary-Safe™ mobile phase safety caps
- Replace every six months for optimum protection.

Cat. No.	Description
CX0001	Replacement Filter Check Valve for Mobile Phase Safety Caps