



analytical

GL38 Caps

2026

Canary + Safe™
Solvent Safety products

GL38
Mobile Phase Bottle Safety Caps



CCGL3801



CCGL3802



CCGL3803



- **Screws directly onto GL38 threaded amber solvent bottles** (NOTE: Most 4L solvent bottles are GL38 threads)
- **No thread adapter needed!**
- **Comes with Air Check Valve and Fittings**

GL38 Mobile Phase Bottle Safety Caps

Cat. No.	Description
CCGL3801	Canary-Safe Cap, GL38, 1 Standard Tubing Port
CCGL3802	Canary-Safe Cap, GL38, 2 Standard Tubing Ports
CCGL3803	Canary-Safe Cap, GL38, 3 Standard Tubing Ports



Replacement Filter Check Valve

- Fits Canary-Safe™ mobile phase safety caps
- Replace every six months for optimum protection.

Cat. No.	Description
CX0001	Replacement Filter Check Valve for Mobile Phase Safety Caps