


Dried Blood Spot Card Punch

DBS Dried Blood Spot Card Punch

The DBS Punch eliminates hand fatigue and blistering from the use of hand-held manual card punches. It is a pneumatic foot-controlled "power" punch that increases speed and efficiency when loading dried blood spot samples into collection plates or tubes.

- Easy to Use - Foot Pedal Operation
- No Hand Fatigue or Blistering
- Reduces Processing Time
- Direct Insertion of Punched Disc into Collection Plate Well or Tube
- New, Better View of Punching Area
- 3mm, 4mm, 5mm and 6mm Punch Sizes Available



DBS Complete Systems

Includes Punch Pod Assembly, Base Stand with Shelf, Foot Pneumatics and Tubing

Cat. No.	Description	Qty
327300	3mm Dried Blood Spot Card Punch, Complete System	Each
327400	4mm Dried Blood Spot Card Punch, Complete System	Each
327500	5mm Dried Blood Spot Card Punch, Complete System	Each
327600	6mm Dried Blood Spot Card Punch, Complete System	Each

Punch Pod Only

Includes Foot Pneumatics and Tubing

Cat. No.	Description	Qty
327301	3mm DBS Punch Pod with Foot Pneumatics and Tubing	Each
327401	4mm DBS Punch Pod with Foot Pneumatics and Tubing	Each
327501	5mm DBS Punch Pod with Foot Pneumatics and Tubing	Each
327601	6mm DBS Punch Pod with Foot Pneumatics and Tubing	Each


Replacement Cutters *(cannot be interchanged between Punch Pod sizes)*

Cat. No.	Description	Qty
327903	Cutter for 3mm DBS Punch	Each
327904	Cutter for 4mm DBS Punch	Each
327905	Cutter for 5mm DBS Punch	Each
327906	Cutter for 6mm DBS Punch	Each

For pricing, see
analytical-sales.com

Replacement Pipettes

Cat. No.	Description	Qty
327888	Glass Pipette for DBS Punch	Each



Operation requires air supply, set between 90-110 psi.