



Fingertight UPLC Fitting

Double Ferrule Stainless Steel Fitting

- Rated to 13,500psi
- Unique Double Ferrule, Fingertight, ZDV
- The Grooved Stainless Steel Collet Reinforces the Polymer Ferrule

Cat. No.	Description	Qty
SDY-201	SS Nut, includes Polymer Ferrule and SS Collet	2
SDY-102	Replacement Polymer Ferrules	5
SDY-103	Replacement Stainless Steel Collets	2

For pricing, see analytical-sales.com