

## Echelon C18: Sub 2 $\mu$ m & 3 $\mu$ m Columns



### ***Echelon C18 provides:***

- Proprietary highest quality C18 & packing method
  - Better retention and resolution of compounds more polar than traditional C18 can manage
  - **Extreme Durability:**
    - pH 2-11
    - Handles 100% aqueous
  - **Proven Reliability:**
    - Adopted by high-throughput screening groups for years
- Featured in 2013 ASMS poster, visit:*  
<http://www.analytical-sales.com/Echelon-Columns.html>

## ***All this now in a UHPLC Column!***

### **<2 $\mu$ m Echelon Columns**

**C18 100Å 1.7 $\mu$ m**

| Cat. No. | Description                                   |
|----------|-----------------------------------------------|
| ADV8106  | Echelon C18, 100Å, 1.7 $\mu$ m, 20mm x 2.1mm  |
| ADV8107  | Echelon C18, 100Å, 1.7 $\mu$ m, 30mm x 2.1mm  |
| ADV8108  | Echelon C18, 100Å, 1.7 $\mu$ m, 50mm x 2.1mm  |
| ADV8117  | Echelon C18, 100Å, 1.7 $\mu$ m, 100mm x 2.1mm |

### **3 $\mu$ m Echelon Columns**

**C18 100Å 3 $\mu$ m**

|         |                                            |
|---------|--------------------------------------------|
| ADV8302 | Echelon C18, 100Å, 3 $\mu$ m, 20mm x 2.1mm |
| ADV8303 | Echelon C18, 100Å, 3 $\mu$ m, 30mm x 2.1mm |
| ADV8304 | Echelon C18, 100Å, 3 $\mu$ m, 50mm x 2.1mm |

*For pricing, see [analytical-sales.com](http://analytical-sales.com)*