

Inserts for 96-Well Collection Plates





96405

96847



961879



* 1mL glass Inserts extend beyond the top surface of the collection plate - use these **round-well** cap mats:

1mL Glass Inserts

Loose & Assembled 1mL Flat Bottom Inserts

Cat. No.	Description	Qty
96405	1mL Flat Bottom Glass Inserts	100
96481	96-Well Plate, 1mL Flat Bottom Glass Inserts	1
96485	96-Well Plate, 1mL Flat Bottom Glass Inserts, 99685 Silicone/PTFE Cap Mat	1
96490	96-Well Plate, 1mL Flat Bottom Glass Inserts, 99690 Pre-slit Silicone/PTFE Cap Mat	1
96489	96-Well Plate, 1mL Flat Bot. Glass Inserts, Pre-slit Opaque Blue Silicone/PTFE Cap Mat	1
96426	96-Well Plate, 1mL Flat Bottom Glass Inserts, Sealing Silicone Cap Mat	1

Loose & Assembled 1mL Tapered Inserts

Cat. No.	Description	Qty
96847	Tapered 1mL Glass Inserts	100
96681	96-Well Plate, Tapered 1mL Glass Inserts	1
99861	96-Well Plate, Tapered 1mL Glass Inserts, Clear Silicone/PTFE Cap Mat	1
99862	96-Well Plate, Tapered 1mL Glass Inserts, Pre-Slit Clear Silicone/PTFE Cap Mat	1
961879	96-Well Plate, Tapered 1mL Glass Inserts, Pre-slit Opaque Blue Silicone/PTFE Cap Mat	1

Loose & Assembled 1mL Flat Bottom Inserts - Silanized

Cat. No.	Description	Qty
96405S	1mL Flat Bottom Tapered Silanized Glass Inserts	100


Loose & Assembled 1mL Tapered Inserts - Silanized

Cat. No.	Description	Qty
96847S	Tapered 1mL Silanized Glass Inserts	100
96681S	96-Well Plate, Tapered 1mL Silanized Glass Inserts	1
99862S	96-Well Plate, Tapered 1mL Silanized Glass Inserts, Pre-Slit Clear Silicone/PTFE Cap Mat, Pre-assembled	1


96-Well Round Cap Mats

Cat. No.	Description	Qty
99685	CLEAR Silicone/PTFE Cap Mat for 1mL Vials	5
99690	CLEAR Silicone/PTFE Pre-Slit Cap Mat for 1mL Vials	5
99890	Blue Opaque Silicone/PTFE Pre-Slit Cap Mat for 1mL Vials	5


Our glass vials are made from 33 expansion borosilicate glass which conforms to ASTM E-438 Type 1-Class A, USP Type 1 and ISO 3585.







99685





99690



99890



* 0.7mL glass Inserts do not reach the top surface of the collection plate - use these **square-well** cap mats:



0.7mL Glass Inserts

Loose & Assembled 0.7mL Flat Bottom Inserts

Cat. No.	Description	Qty
96305	0.7mL Flat Bottom Glass Inserts	100
96382	96-Well Plate, 0.7mL Flat Bottom Glass Inserts	1
96326	96-Well Plate, 0.7mL Flat Bottom Glass Inserts, 2mL Polypropylene Pierceable Cap Mat with "X" on Surface	1
96380	96-Well Plate, 0.7mL Flat Bottom Glass Inserts, Pre-Slit Silicone/PTFE Cap Mat, Pre-Assembled	1
96385	96-Well Plate, 0.7mL Flat Bottom Glass Inserts, CLEAR Ultra Thin Silicone/PTFE Cap Mat, Pre-Assembled	1

Loose & Assembled 0.7mL Tapered Inserts

Cat. No.	Description	Qty
96837	Tapered 0.7mL Glass Inserts	100
96687	96-Well Plate, Tapered 0.7mL Glass Inserts	1
96627	96-Well Plate, Tapered 0.7mL Glass Inserts, Polypropylene Cap Mat	1
96685	96-Well Plate with Tapered 0.7mL Glass Inserts and Purple Ultra Thin Pre-Slit Square Cap Mat, Pre-assembled	1

Loose & Assembled 0.7mL Flat Bottom Inserts - Silanized

Cat. No.	Description	Qty
96305S	0.7mL Flat Bottom Silanized Glass Inserts	100
96326S	0.7mL Flat Bottom Silanized Glass Inserts with Cap Mat & Plate	1

Loose & Assembled 0.7mL Tapered Inserts - Silanized

Cat. No.	Description	Qty
96837S	Tapered 0.7mL Silanized Glass Inserts	100
96687S	96-Well Plate, Tapered 0.7mL Silanized Glass Inserts	1
96627S	96-Well Plate with Silanized Tapered 0.7mL Glass Inserts and Pierceable Polypropylene Cap Mat, Pre-assembled	1

96-Well Square Cap Mats

Cat. No.	Description	Qty
2726	2mL Autosampler Compatible Square Polypropylene Cap Mat with 'X' on Surface	20
964008	CLEAR Ultra Thin Silicone/PTFE Square Cap Mat	5
964085	CLEAR Pre-Slit , Ultra Thin Silicone/PTFE Square Cap Mat	5
965008	Purple Ultra Thin Silicone/PTFE Square Mighty Mat	5
965085	Purple Pre-Slit , Ultra Thin Silicone/PTFE Square Mighty Mat	5

For pricing, see analytical-sales.com