




Lumidox Controller

LUM96B
LED Array (sold separately)


Lumidox System

LED Photon Generator for Photoredox Applications


Intelligent Autosensing LED controller for 24 well and 96 well LED arrays



LUM96G



LUM96W



LUM96BGW

- Repeatability of milliamp selection for photon intensity
- Non-magnetic LED base
- Robust Controller Shell
- PTFE coated cabling
- Chemically resistant LED cover
- 30mA LED Maximum Output
- 470nm, 527nm, White LED's
- Available in **Blue**, **Green**, White and Alternating

Catalog No.	Description
LUM24B	Lumidox 24-Well Blue LED Array
LUM24G	Lumidox 24-Well Green LED Array
LUM96G	Lumidox 96-Well Green LED Array
LUM96W	Lumidox 96-Well White LED Array
LUM96BGW	Lumidox 96-Well, Alternating Blue , Green and White LED Array
LUMCON	Lumidox LED Controller (<i>LED arrays sold separately</i>)
UVGoggles	UV Safety Goggles

For pricing, see analytical-sales.com




Lumidox® UV LED System



Lumidox UV Controller



LUM96-365
UV LED Array (sold separately)



Lumidox UV System

UV LED Photon Generator for Photoredox Applications

Intelligent Autosensing LED controller for 96 well UV LED arrays

- Repeatability of milliamp selection for photon intensity
- Non-magnetic LED base
- Robust Controller Shell
- PTFE coated cabling
- Chemically resistant LED cover
- 20mA LED Maximum Output
- 365nm and 415nm UV LED's
- Includes 1 pair of UV protective goggles
- *Use Lumidox UV Controller with UV Arrays Only!*

Catalog No.	Description	Qty
LUMCON-UV	Lumidox UV LED Controller for Ultraviolet 365nm and Violet 415nm Arrays (<i>arrays sold separately</i>)	Each
LUM96-365	LumidoxUV 96-Well, 365nm, Ultraviolet LED Array	Each
LUM96-415	LumidoxUV 96-Well, 415nm, Violet LED Array	Each
UVGoggles	Additional UV Safety Goggles	Each

For pricing, see analytical-sales.com



The Lumidox UV System comes with one pair of UV protective goggles. Additional goggles can be purchased.