



Filterless (or Fritless) Plate for Solid Phase Extraction (SPE)

- No filters, glues, polymers, or matrices
- No sample loss or contamination risk on filter
- 300 μ L well volume for 96-well plates
- SPE bed sizes for either 2-3 μ L or 10 μ L of sample

| Cat. No. | Description |
|----------|-------------|
|----------|-------------|

| | |
|-------|---|
| 96311 | Frit-less Micro-elution (2-3 μ L) 96 well SPE plate, Silica WCX |
|-------|---|

| | |
|-------|--|
| 96312 | Frit-less Micro-elution (10 μ L) 96 well SPE plate, Silica WCX |
|-------|--|

| | |
|-------|---|
| 96321 | Frit-less Micro-elution (2-3 μ L) 96 well SPE plate, Silica WAX |
|-------|---|

| | |
|-------|--|
| 96322 | Frit-less Micro-elution (10 μ L) 96 well SPE plate, Silica WAX |
|-------|--|

For pricing, see analytical-sales.com

Also available in 384 well versions (30 μ L well volume).

Many other SPE chromatographic media available:

Silica, C-18, C-8, C-4, Carbon, Zirconium dioxide, Titanium dioxide, Protein A, Protein G, Gel Filtrations, etc...