

MonoSLEEVE™



MonoSLEEVE™ Single Column Heater Controller

- Self-calibrating
- Plug & Play smart heater technology
- Temperature range from 25°C to 85°C
- Controller comes with Power Supply, USB Cable, GUI Software and Manual.



MonoSLEEVE is the next generation of the original HotSleeve controller. With it's advanced "Plug & Play" technology, MonoSleeve can accommodate single **HotSleeve™+*** or **AgileSleeve+®*** column heaters of different sizes without having to recalibrate to each individual sleeve. This allows you to switch out different heaters on the fly.

** HotSleeve+ is the "smart" version of the original HotSleeve. AgileSleeve+ is the "smart" version of AgileSleeve.*

Cat. No.	Description
MonoSLEEVE	Single Column Heater Controller. Includes: MonoSLEEVE Controller, Power Supply, USB Cable, GUI, and Manual. <i>Heater not included.</i>

Note: MonoSLEEVE will only accept HotSleeve+ and AgileSleeve+ heaters with "Plug & Play" technology.
See information about these heaters on back side.

HotSleeve™ + and AgileSleeve+®

All heaters made by Analytical Sales & Services are calibrated to NIST traceable standards.

HotSleeve+ and AgileSleeve+ are self-calibrating “smart” heaters that utilize “Plug & Play” technology. This technology eliminates error generated when switching heaters from one system to another.

HotSleeve™+

Analytical Size Column Heater

- Self-calibrating
- Temperature range from 5°C above ambient to 85°C

Cat. No.	Description
----------	-------------


HSI-5L	HotSleeve+ Heater for 5cm Columns
--------	-----------------------------------

HSI-10L	HotSleeve+ Heater for 10cm Columns
---------	------------------------------------

HSI-15L	HotSleeve+ Heater for 15cm Columns
---------	------------------------------------

HSI-25L	HotSleeve+ Heater for 25cm Columns
---------	------------------------------------

HSI-30L	HotSleeve+ Heater for 30cm Columns
---------	------------------------------------



AgileSleeve+®

Flexible Capillary Heater

- Accurate from 5°C above ambient to 100°C
- Heat lines up to 30cm long, with outside diameters of 1/32" or 1/16"
- Great for small diameter analytical columns, nanospray, injector or pre-column mobile phase lines

Cat. No.	Description
----------	-------------

ASI1032H	AgileSleeve+ 10cm x 1/32" ID Self-calibrating Heater
----------	--

ASI1016H	AgileSleeve+ 10cm x 1/16" ID Self-calibrating Heater
----------	--

ASI1532H	AgileSleeve+ 15cm x 1/32" ID Self-calibrating Heater
----------	--


ASI1516H	AgileSleeve+ 15cm x 1/16" ID Self-calibrating Heater
----------	--

ASI2532H	AgileSleeve+ 25cm x 1/32" ID Self-calibrating Heater
----------	--

ASI2516H	AgileSleeve+ 25cm x 1/16" ID Self-calibrating Heater
----------	--

ASI3032H	AgileSleeve+ 30cm x 1/32" ID Self-calibrating Heater
----------	--

ASI3016H	AgileSleeve+ 30cm x 1/16" ID Self-calibrating Heater
----------	--



All Analytical Sales and Services heater products carry the value recognized CE mark. Each heater meets the demanding electromagnetic emission standards of the new European Union Directives, United States standards, and Canadian standards.

For pricing, see analytical-sales.com



179 Rt 206, Flanders, NJ 07836

Phone: 973-616-0700 • Fax: 973-616-0133

info@analytical-sales.com • analytical-sales.com

©2023 Analytical Sales and Services, Inc.

Doc. Ref.: MonoSLEEVE rev1d