




Multi-Inlet Waste Manifold



Keep your solvent waste lines neat and organized with our new Multi-Inlet Waste System Manifold!


- Made of durable chemical-resistant polypropylene
- Fits bench up to 3" thick
- Includes fittings and barbed adapters
- Other configurations available - call for details



For 1/16" or 2.5mm tubing, use these fittings (sold separately)





| Cat. No. | Description | Qty |
|----------|---|-----|
| 59249 | PEEK™ One Piece Fitting (1/4-28) for 1/16" OD Tubing | Ea |
| 59305X | Green Delrin Flangeless Male Nut, use with Ferrule (59353X) | 10 |
| 59353X | ETFE Ferrule for 2.5mm OD Tubing | 10 |



Plugs

| Cat. No. | Description | Qty |
|----------|--|-----|
| CF5007 | Plugs For 1/8" or 1/16" Ports, White PFA | 10 |
| CF5009 | Plugs for 6.4mm ID Barbed Connector Port, PP | 5 |



Multi-Inlet Waste Manifold

| Cat. No. | Description |
|----------|---|
| MW0802 | Canary-Safe Multi-Inlet Waste Manifold with 8 Standard Tubing Ports for 1/16" or 1/8" OD Tubing, 2 Ports for Barbed Adapter, and 1 Large Outlet. Includes 1/8" Fittings and Barbed Adapters |

Other configurations available, see analytical-sales.com

For pricing, see analytical-sales.com

The Canary-Safe™ **Multi-inlet Waste Manifold** remediates these waste solvent situations:



Exhaust Filter Adapters for JustRite® Waste Containers

| Cat. No. | Description | Qty |
|----------|---|------|
| CJX01 | Exhaust Filter Fitting/ <i>Long</i> Adapter for SS Port on JustRite Waste Container | Each |
| CJX02 | Exhaust Filter Fitting/ <i>Short</i> Adapter for Gray PEEK Port on JustRite Waste Container | Each |



- Timer/Change Indicator**
- 1. Push button to activate
 - 2. Replace when service life expires

Exhaust Filters with Timer/Change Indicator

Canary-Safe exhaust filters contain an activated carbon, multi-component granulate material which absorbs 99% of all vapor and particulate materials evaporating from your solvents. Filters come with a **Timer/Change Indicator** at no additional cost.

Canary-Safe™ Vapor Exhaust Filters for Safety Waste Caps

| Cat. No. | Size | Description | Replacement | Qty |
|----------|------|--|-----------------|-----|
| CX3000 | S | Exhaust Filter, for containers < 5L | Every 3 Months | Ea |
| CX6000 | M | Exhaust Filter, for containers < 20L | Every 6 Months | Ea |
| CX12000 | L | Exhaust Filter, for containers > 20L | Every 12 Months | Ea |
| CX6010 | M | Exhaust Filter with H2 SO4 (sulfuric acid) for Volatile Bases | Every 6 Months | Ea |
| CX6020 | M | Exhaust Filter with KOH (potassium hydroxide) for Volatile Acids | Every 6 Months | Ea |

