



MultiSLEEVE™

Multiple Column Heater Controller

- Self-calibrating
- Plug & Play smart heater technology
- Controls up to 4 column heaters simultaneously
- 4 independent channels
- Used in multiplexing applications
- Temperature range from 5°C above ambient to 85°C
- Controller comes with cables, manual and USB GUI software.

The MultiSLEEVE controller is used for applications which require multiple heating zones or multiple LC system heating. Used in conjunction with HotSleeve+™ and/or AgileSleeve+® heaters (see reverse), it maintains four independent channels simultaneously and accurately while reducing overall equipment costs.

Note: MultiSLEEVE will only accept HotSleeve+ and AgileSleeve+ heaters with “Plug & Play” technology. See information about these heaters on back side.



All Analytical Sales and Services heater products carry the value recognized CE mark. Each heater meets the demanding electromagnetic emission standards of the new European Union Directives, United States standards, and Canadian standards.

| Cat. No. | Description |
|-------------|-----------------------------------|
| MultiSLEEVE | Multiple Column Heater Controller |

For pricing, see analytical-sales.com

HotSleeve+™ and AgileSleeve+®

All heaters made by Analytical Sales & Services are calibrated to NIST traceable standards.

HotSleeve+ and AgileSleeve+ are self-calibrating “smart” heaters that utilize “Plug & Play” technology. This technology reduces all of the error generated when switching heaters from one system to another.

HotSleeve+™

Analytical Size Column Heater

- Self-calibrating
- Temperature range from 5°C above ambient to 85°C

| Cat. No. | Description |
|----------|---|
| HSI-5L | HotSleeve+ Self-calibrating Heater for 5cm Columns |
| HSI-10L | HotSleeve+ Self-calibrating Heater for 10cm Columns |
| HSI-15L | HotSleeve+ Self-calibrating Heater for 15cm Columns |
| HSI-25L | HotSleeve+ Self-calibrating Heater for 25cm Columns |
| HSI-30L | HotSleeve+ Self-calibrating Heater for 30cm Columns |

Call to inquire about Prep, Semi-prep and Super size HotSleeve Plus



AgileSleeve+™

Flexible Capillary Heater

- Accurate from 30°C to 100°C
- Heat lines up to 30 cm long, with outside diameters of 1/32" or 1/16"
- Great for small diameter analytical columns, nanospray, injector or pre-column mobile phase lines

| Cat. No. | Description |
|----------|---|
| ASI1032H | AgileSLEEVE+ 10cm x 1/32" ID Self-calibrating Heater |
| ASI1016H | AgileSLEEVE+ 10cm x 1/16" ID Self-calibrating Heater |
| ASI1532H | AgileSLEEVE+, 15cm x 1/32" ID Self-calibrating Heater |
| ASI1516H | AgileSLEEVE+ 15cm x 1/16" ID Self-calibrating Heater |
| ASI2532H | AgileSLEEVE+ 25cm x 1/32" ID Self-calibrating Heater |
| ASI2516H | AgileSLEEVE+ 25cm x 1/16" ID Self-calibrating Heater |
| ASI3032H | AgileSLEEVE+ 30cm x 1/32" ID Self-calibrating Heater |
| ASI3016H | AgileSLEEVE+ 30cm x 1/16" ID Self-calibrating Heater |

For pricing, see analytical-sales.com



All Analytical Sales and Services heater products carry the value recognized CE mark. Each heater meets the demanding electromagnetic emission standards of the new European Union Directives, United States standards, and Canadian standards.