



## MultiSLEEVE™

### Multiple Column Heater Controller

- Self-calibrating
- Plug & Play smart heater technology
- Controls up to 4 column heaters simultaneously
- 4 independent channels
- Used in multiplexing applications
- Temperature range from 5°C above ambient to 85°C
- Controller comes with cables, manual and USB GUI software.

The MultiSLEEVE controller is used for applications which require multiple heating zones or multiple LC system heating. Used in conjunction with HotSleeve+™ and/or AgileSleeve+® heaters (see *reverse*), it maintains four independent channels simultaneously and accurately while reducing overall equipment costs.

**Note: MultiSLEEVE will only accept HotSleeve+ and AgileSleeve+ heaters with “Plug & Play” technology. See information about these heaters on back side.**



All Analytical Sales and Services heater products carry the value recognized CE mark. Each heater meets the demanding electromagnetic emission standards of the new European Union Directives, United States standards, and Canadian standards.

Cat. No.	Description
MultiSLEEVE	Multiple Column Heater Controller

For pricing, see [analytical-sales.com](http://analytical-sales.com)

# HotSleeve+™ and AgileSleeve+®

All heaters made by Analytical Sales & Services are calibrated to NIST traceable standards.

HotSleeve+ and AgileSleeve+ are self-calibrating “smart” heaters that utilize “Plug & Play” technology. This technology reduces all of the error generated when switching heaters from one system to another.


## HotSleeve+™

### Analytical Size Column Heater

- Self-calibrating
- Temperature range from 5°C above ambient to 85°C

Cat. No.	Description
HSI-5L	HotSleeve+ Self-calibrating Heater for 5cm Columns
HSI-10L	HotSleeve+ Self-calibrating Heater for 10cm Columns
HSI-15L	HotSleeve+ Self-calibrating Heater for 15cm Columns
HSI-25L	HotSleeve+ Self-calibrating Heater for 25cm Columns
HSI-30L	HotSleeve+ Self-calibrating Heater for 30cm Columns

Call to inquire about Prep, Semi-prep and Super size HotSleeve Plus



## AgileSleeve+™

### Flexible Capillary Heater

- Accurate from 30°C to 100°C
- Heat lines up to 30 cm long, with outside diameters of 1/32” or 1/16”
- Great for small diameter analytical columns, nanospray, injector or pre-column mobile phase lines

Cat. No.	Description
ASI1032H	AgileSLEEVE+ 10cm x 1/32” ID Self-calibrating Heater
ASI1016H	AgileSLEEVE+ 10cm x 1/16” ID Self-calibrating Heater
ASI1532H	AgileSLEEVE+, 15cm x 1/32” ID Self-calibrating Heater
ASI1516H	AgileSLEEVE+ 15cm x 1/16” ID Self-calibrating Heater
ASI2532H	AgileSLEEVE+ 25cm x 1/32” ID Self-calibrating Heater
ASI2516H	AgileSLEEVE+ 25cm x 1/16” ID Self-calibrating Heater
ASI3032H	AgileSLEEVE+ 30cm x 1/32” ID Self-calibrating Heater
ASI3016H	AgileSLEEVE+ 30cm x 1/16” ID Self-calibrating Heater

For pricing, see [analytical-sales.com](http://analytical-sales.com)



All Analytical Sales and Services heater products carry the value recognized CE mark. Each heater meets the demanding electromagnetic emission standards of the new European Union Directives, United States standards, and Canadian standards.