

A Selection of Products for Your Workflow - From Sample to Detector



24108



24727



24777

Collection Plates & Cap Mats

24-Well Square Reservoir Plates and Cap Mat

Cat. No.	Format	Description	Material	Bottom	Volume	Qty
24108	24-Well	Square Well Reservoir Plate	Clear PP	Pyramid	10mL	25
24108GC	24-Well	Square Well Reservoir Plate w/ SiliGuard	Clear PP	Pyramid	10mL	25
24727	24-Well	Square Well Reservoir Plate	Clear PP	Round	10mL	25
24727GC	24-Well	Square Well Reservoir Plate w/ SiliGuard	Clear PP	Round	10mL	25
24777	24-Well	CLEAR Silicone Cap Mat				5



59623-23

Advantage™ 2mL 96-Well Collection Plates

Square Well with V-Shaped Bottom


Cat. No.	Description	Material	Volume	Qty
59623-23	Square Top, V-Bottom Collection Plate	PP	2mL	10
59623-23GC	Square Top, V-Bottom w/ SiliGuard	PP	2mL	10

Resealable Cap Mats

- Compatible with all current autosampler needles
- No coring or cross talk
- Fits all standard 96-position blocks
- Reseals after multiple injections
- Pre-slit or non-pre-slit

Cat. No.	Description	Color	Qty
965085	Pre-Slit , Ultra Thin Square Well, Soft Silicone/PTFE	Purple	5

965085




Other items, options and/or sizes are available - please call for details.


Sealing Film + (Patented) Pattern Adhesive PTFE (Teflon®) Sealing Film

- Adhesive-free target areas
- Firm seal to microplate -80°C to +120°C
- Prevents cross-contamination and evaporation

- Outstanding performance with Autosamplers
 - Only 0.5 lb force needed to pierce
 - No coring of material or clogging of needle
- Optically clear for inspection and diagnosis

Cat. No.	Description	Qty
961201	Advantage™ Teflon (PTFE) Adhesive Film, Square 96-Well Pattern, Ultra-thin	100





Vacuum Manifold Filtration System for SPE Sample Preparation


- Adjustable vacuum control
- Compatible with all standard filter bottom microplates
- Collect into any storage plate
- Compatible with robotic handling

Cat. No. Description

96844	Vacuum Manifold Filtration System for SPE Sample Filtration
-------	---

24-well Filter Plates

Cat. No.	Description	Qty
24152	24 Well Filter Plates, PE 25µm, 15mL	15




Canary Safe™ Solvent Safety Products

Mobile Phase Bottle Safety Caps

Cat. No. Description

CC1001	Canary-Safe™ Cap 1, GL45, 1 Standard Port for 1/8" OD Tubing
CC1002	Canary-Safe™ Cap 2, GL45, 2 Standard Ports for 1/8" OD Tubing

UHPLC Finger Tight Assemblies



Ultra High Pressure HPLC Fitting Assemblies are the perfect column connectors for HPLC because they totally eliminate any possibility of Dead Volume in the critical inlet and outlet areas, while holding pressure to **27,000 psi**.

Single UHPLC Assembly

Cat. No.	Dimensions
61650	50cm x 0.010"ID
61250	50cm x 0.020"ID

Double UHPLC Assembly

Cat. No.	Dimensions
62650	50cm x 0.010"ID
62250	50cm x 0.020"ID

All PEEK Tubing ID & Lengths also Available



Lite-Touch™

- Rated to **5,000 psi**



For 1/16" OD Tubing

Cat. No.	Description	Qty
68561-1	PEEK™ Ferrule, SS Ring	5
68562-1	PEEK™ Nut, 10-32	5

For 1/8" OD Tubing

Cat. No.	Description	Qty
68563-1	PEEK™ Ferrule, SS Ring	5
68564-1	PEEK™ Nut, 1/4-28	5



One Piece Ferrule System

- Optional replacement for Lite-Touch PEEK ferrules
- Rated Up to **13,500 psi** when used with Lite-Touch

Cat. No.	Description	Qty
68566-1	One Piece Ferrule System	5

For pricing, see analytical-sales.com

This is a quick guide of our top selling products.

Prep Inline Filter

- Replace o-ring periodically
- Use supplied fittings to make proper connections

Cat. No.	Description	Qty
49230	Prep Inline Filter Assembly w/ 2µm Frit	1
49230DR	Prep Inline Filter Assembly w/ 2µm Frit, Larger Flow Ports	1
49231X	Replacement Inline Prep Filter, 2µm	10
49233	Replacement O-Ring for Prep Inline Filter System	1
49501	SS Frit, 0.839" x 0.125" , 10µm	20



MonoSLEEVE™

SINGLE Column Heater Controller

- Self-calibrating
- Plug & Play smart heater technology
- Temperature range from 25°C to 85°C
- Controller comes with Power Supply, USB Cable, GUI Software and Manual.

MonoSLEEVE is the next generation of the original HotSLEEVE controller. With it's advanced "Plug & Play" technology, MonoSleeve can accommodate single **HotSleeve™+*** or **AgileSleeve+®*** column heaters of different sizes without having to recalibrate to each individual sleeve. This allows you to switch out different heaters on the fly.

** HotSleeve Plus is the new "smart" version of the original HotSleeve.
AgileSleeve Plus is the new "smart" version of AgileSleeve.*

Cat. No.	Description
MonoSLEEVE	Single Column Heater Controller. <i>Heater not included.</i>

HotSleeve+™

Prep & Semi-Prep Size Column Heaters

Cat. No.	Description
HSI-30L	HotSleeve+ Heater for 30cm Columns
HSI-1012L	HotSleeve+ Heater - 1" ID x 12" Length
HSI-1512L	HotSleeve+ Heater - 1.5" ID x 12" Length
HSI-2012L	HotSleeve+ Heater - 2" ID x 12" Length

0.80", 1.75" & 3.5" ID are also available.




MultiSLEEVE™

MULTIPLE Column Heater Controller

- Self-calibrating
- Plug & Play smart heater technology
- Controls up to 4 column heaters simultaneously
- 4 independent channels
- Used in multiplexing applications
- Temperature range from 25°C to 85°C
- Controller comes with cables, manual and USB GUI software.

The MultiSLEEVE controller is used for applications which require multiple heating zones or multiple LC system heating. Used in conjunction with **HotSleeve™+*** and/or **AgileSleeve+®*** heaters, it maintains four independent channels simultaneously and accurately while reducing overall equipment costs.

Cat. No.	Description
MultiSLEEVE	Multiple Column Heater Controller. <i>Heaters not included.</i>



Other items, options and/or sizes are available - please call for details.

