




PTFE Reservoirs



NEW! PTFE Reservoirs

- Made from a solid PTFE block
- Precision machined wells, SBS compliant
- Ideal for aggressive organics
- Enables your automation's full pipetting potential/efficiency
- Reusable when used with clean solvents (no analytes, etc)
- Choice of three well configurations (two w/5mL and one w/10mL wells)

Cat. No.	Description	Qty
40805	PTFE Reservoir with Conical Wells, 4 Rows, 2 Columns, 5 mL / Well	Each
40810	PTFE Reservoir with Conical Wells, 4 Rows, 2 Columns, 10 mL / Well	Each
41605	PTFE Reservoir with Conical Wells, 8 Rows, 2 Columns, 5 mL / Well	Each

For pricing, see analytical-sales.com