





SpeedyBead Dropper



NEW! SpeedyBead™ 96-Well Bead Dropper

- Holds up to 1,000 beads
- Fill up to 10 collection plates before reloading
- Drops one agitator bead into each of the 96 wells simultaneously
- Increases efficiency - reduces prep time
- 100% accuracy - no missed wells, no double-up
- For beads up to 5/32" (0.15625")

See it in Action!

Watch the video at:

<http://www.analytical-sales.com/SpeedyBead-Video.html>

Cat. No.	Description
96T051	Speedy-Bead 96-Well Bead Dropper

For pricing, see analytical-sales.com