




analytical


# SpeedyBead Dropper

2024



## **NEW!** SpeedyBead™ 96-Well Bead Dropper

- Holds up to 1,000 beads
- Fill up to 10 collection plates before reloading
- Drops one agitator bead into each of the 96 wells simultaneously
- Increases efficiency - reduces prep time
- 100% accuracy - no missed wells, no double-up
- For beads up to 5/32" (0.15625")



### **See it in Action!**

Watch the video at:

<http://www.analytical-sales.com/SpeedyBead-Video.html>

Cat. No.	Description
96T051	Speedy-Bead 96-Well Bead Dropper

For pricing, see [analytical-sales.com](http://analytical-sales.com)

179 Rt 206, Flanders, NJ 07836

Phone: 973-616-0700 • Fax: 973-616-0133

[customerservice@analytical-sales.com](mailto:customerservice@analytical-sales.com) • [www.analytical-sales.com](http://www.analytical-sales.com)

©2024 Analytical Sales and Services, Inc.

Doc. Ref.: SpeedyBead Dropper rev1e