



## Stainless Steel Top Cover for Reaction Blocks

- 96 Well Reactor Lid with Enlarged (4.3mm) Holes
- 316 Stainless Steel
- Made for Automated Sampling Instruments

**Made in the USA**

Cat. No.	Description	Thickness	Qty
96985	SS 96-Well Reaction Block Cover with Enlarged (4.3mm) Holes	6mm	Each
96986	SS 96-Well Reaction Block Cover with Enlarged (4.3mm) Holes	3mm	Each

*For pricing, see [analytical-sales.com](http://analytical-sales.com)*