




analytical

Stainless Steel Reaction Block Cover

2026



Stainless Steel Top Cover for Reaction Blocks

- 96 Well Reactor Lid with Enlarged (4.3mm) Holes
- 316 Stainless Steel
- Made for Automated Sampling Instruments

Made in the USA

Cat. No.	Description	Thickness	Qty
96985	SS 96-Well Reaction Block Cover with Enlarged (4.3mm) Holes	6mm	Each
96986	SS 96-Well Reaction Block Cover with Enlarged (4.3mm) Holes	3mm	Each

For pricing, see analytical-sales.com