





Thermal Adapter Plates



967720TAP



967720TAP



↓
Place on top of heat block or chiller
(not included)

NEW! Thermal Adapter Plates for Polypropylene Collection Plates

- For improved temperature transfer to samples

The Thermal Adapter Plate conforms to the bottom shape of the collection plate wells to ensure uniform thermal transfer for heating and cooling when used in conjunction with a heating block or chiller (not included). The plate features a SLAS footprint and easy to identify purple anodized coating.

Analytical Sales & Services can make Thermal Adapter Plates to accommodate select polypropylene collection plates (minimum quantities may apply). Please call 973-616-0700 for more information.

Cat. No.	Description	Qty
967720TAP	Thermal Adapter Plate for 967720 PP Collection Plate	Each
59623-23TAP	Thermal Adapter Plate for 59623-23 PP Collection Plate	Each

For pricing, see analytical-sales.com