



1mL TrueTaper® Collection Plates

1mL TrueTaper

1mL, 96-Well Polypropylene Collection Plate


- Round well with unique **tapered bottom** for maximum sample recovery
- 1.0mL total volume capacity
- 50µL tapered bottom volume

Cat. No.	Description	Volume	UOM
968810	1mL TrueTaper 96-Well Collection Plate	50µL bottom/1.0mL total	25
968810GC	1mL TrueTaper w/ SiliGuard	50µL bottom/ 1.0mL total	25


Use with:

Purple Silicone Cap Mats for 1mL TrueTaper Collection Plates

- Unique O-ringed Plugs
- Holds Seal with Vigorous Vortexing
- No Crosstalk
- Autosampler Compatible



Cat. No.	Description	UOM
960181	Purple Silicone/PTFE Cap Mat for 1mL TrueTaper Collection Plate	5
960185	Purple Pre-Slit Silicone/PTFE Cap Mat for 1mL TrueTaper Plate	5






50µL well bottom volume

Or use with:

96-Well Round Cap Mats

Cat.No.	Description	UOM
99685	Advantage CLEAR Silicone/PTFE Cap Mat for 1mL TrueTaper Plate and 1mLTapered Vials	5
99690	Advantage CLEAR Pre-Slit Silicone/PTFE Cap Mat for 1mL TrueTaper Plate and 1mL Tapered Vials	5
99890	Advantage Blue Opaque Pre-Slit Silicone/PTFE Cap Mat for 1mL TrueTaper Plate and 1mLTapered Vials	5



Or use with:

Silicone Film

Pattern Adhesive Resealable Silicone Sealing Film (patented)


Cat. No.	Description	UOM
96104	Round 96-Well Adhesive Silicone Sealing Film	5

Teflon Films

Pattern Adhesive PTFE (Teflon®) Sealing Film (patented)

Cat. No.	Description	UOM
961801	Advantage™ Teflon (PTFE) Adhesive Film, Round 96-Well Pattern, Ultra-thin	100
961802	Advantage™ Pre-scored Teflon (PTFE) Adhesive Film, Round 96-Well Pattern, Ultra-thin	100

For pricing, see analytical-sales.com



Color may vary

2mL TrueTaper® Collection Plates

2mL TrueTaper

2mL, 96-Well Collection Plate - Square Well



- Square well with unique **tapered bottom** well design for maximum sample recovery
- 2.0mL total volume capacity
- 100µL tapered bottom volume



Cat. No.	Description	Material	Volume	UOM
968820	TrueTaper 96-Well Collection Plate	PP	100µL bottom/2.0mL total	25
968820GC	TrueTaper w/ SiliGuard	PP	100µL bottom/ 2.0mL total	25

Black TrueTaper for Light Sensitive Compounds

Cat. No.	Description	Material	Volume	UOM
968820BL	TrueTaper 96-Well Collection Plate, Black	PP	100µL bottom/2.0mL total	25



963724

Standard Square Cap Mats

- Thin polypropylene – less than 1/16" at the "X"
- Prevents cross contamination between wells • Economical

Cat. No.	Description	Color	UOM
2726	2mL Autosampler Compatible Square Cap Mat, Pierceable 'X' on Surface	CLEAR	20
3726	2mL Square Cap Mat, Low Temp. (-80°), For Long Term Storage	CLEAR	24




964085

Resealable Cap Mats

- Compatible with all current autosampler needles
- No coring or cross talk • Fits all standard 96-position blocks
- Reseals after multiple injections • Pre-slit or non-pre-slit

Cat. No.	Description	Color	UOM
963725	Square Well, Silicone/PTFE	Green	5
963724	Pre-Slit , Square Well, Silicone/PTFE	Green	5
964008	Ultra Thin Square Well, Soft Silicone/PTFE	CLEAR	5
964085	Pre-Slit , Ultra Thin Square Well, Soft Silicone/PTFE	CLEAR	5
965008	Ultra Thin Square Well, Soft Silicone/PTFE	Purple	5
965085	Pre-Slit , Ultra Thin Square Well, Soft Silicone/PTFE	Purple	5
29660	Square Well, Soft Silicone/PTFE, good with aggressive solvents and extreme shaking	Gray	5
29665	Pre-Slit , Square Well, Soft Silicone/PTFE, good w/ aggressive solvents	Gray	5



965085




29660

Thicker Resealable Silicone, PTFE Coated Cap Mats

- Recommended for Shimadzu autosamplers
- Better suited for 2° to 8°C autosampler temperatures

Cat. No.	Description	Color	UOM
968802	Square Well, Soft Silicone/PTFE	Blue	5
968801	Pre-Slit , Square Well, Soft Silicone/PTFE	Blue	5

968801



961201

Adhesive Teflon Sealing Film

with Adhesive-free Well Patterns for Collection Plates

Cat. No.	Description	UOM
961201	Advantage™ Teflon (PTFE) Adhesive Film, Square 96-Well Pattern, Ultra-thin	100
961202	Advantage™ Pre-scored Teflon (PTFE) Adhesive Film, Square 96-Well Pattern, Ultra-thin	100