


O-ringed Plugs

Purple Silicone Cap Mats for 1mL TrueTaper Collection Plates

- Unique O-ringed Plugs
- Holds Seal with Vigorous Vortexing
- No Crosstalk
- Autosampler Compatible

Cat. No.	Description	Qty
960181	Purple Silicone/PTFE Cap Mat for 1mL TrueTaper Collection Plate	5
960185	Purple <i>Pre-Slit</i> Silicone/PTFE Cap Mat for 1mL TrueTaper Collection Plate	5

1mL TrueTaper
Collection Plate
(968810)Purple Silicone Cap Mat
(960185)

50µL well bottom volume

1mL TrueTaper Collection Plates

Cat. No.	Description	Volume	Qty
968810	1mL TrueTaper™ 96-Well Collection Plate	50µL bottom/1.0mL total	25
968810GC	1mL TrueTaper™ with SiliGuard Low-bind Coating	50µL bottom/ 1.0mL total	25