



## X-Type Syringes for CTC PAL<sup>®</sup> LC Autosamplers

### • Deactivated Needle and Barrel

- Eliminates contact between samples and metal components
- Prevents sample adsorption on inner surface of the barrel
- Extremely smooth surface for low friction and increased plunger lifetime

### • Inert PTFE-Polymer Plunger Tip

- High mechanical stability for increased syringe lifetime
- Highly inert plunger tip material compatible with organics and high and low pH solutions

### • Unique Needle Design

- Needle attached directly to the barrel
- No glue or cement in fluid path
- Improved design for limited carryover

#### Fixed Needle (FN) Syringe (point style 3)

Cat. No.	Volume	Model	Gauge	CTC Part No.
42044754	25 $\mu$ L	1702	22s	SYRX G25-22S-3
42044004	100 $\mu$ L	1710	22	SYRX G100-22-3
42044524	100 $\mu$ L	1710	22s	SYRX G100-22S-3

For pricing, see [analytical-sales.com](http://analytical-sales.com)