





- NOTES:
- 1. MATERIAL: WHITE PTFE (TEFLON)
 - 2. PLATE CONSISTS OF 16 WELLS
 - 3. VOLUME OF EACH WELL NOMINALLY 5mL




SECTION A-A



DETAIL B
SCALE 4 : 1

REVISION HISTORY			
REV	DESCRIPTION	DATE	APPROVED
1A	INITIAL RELEASE	7/1/2019	ENGINEERING



P

DRAWN Matt Dabrowski email engineering@analytical-sales.com	7/1/2019	ANALYTICAL SALES & SERVICES, INC © 2019		
		TITLE		
		41605_16well_5mL_rev1a_P		
	SIZE D			REV 1A
	SCALE 3 : 2		SHEET 1 OF 1	