



NOTES:

1. MATERIAL: ANODIZED ALUMINUM PLATE WITH KAPTON LINER
2. UPPER LEFT CORNER NOTCHED FOR QUICK A1 IDENTIFICATION
3. APPROX. VOLUME: 400 μ L/WELL
4. PLATE WILL BE ENGRAVED WITH ALPHANUMERICS

REVISION HISTORY			
REV	DESCRIPTION	DATE	APPROVED
1A	INITIAL RELEASE	4/22/2026	ENGINEERING



DRAWN MD	4/22/2026	ANALYTICAL SALES AND SERVICES, INC © 2026	
email engineering@analytical-sales.com			
TITLE		482540_rev1a_PUBLIC	
SIZE B			REV 1A
SCALE	1 : 1	SHEET 1 OF 1	

P