

NOT TO SCALE

127.76 ±0.25

99.00 ±0.25

14.38 ±0.25

9.00

Y

27.00 ±0.10

0.50

SECTION Y-Y

ø7.00 ±0.10
(x96)

4.00 (x8)

25.00

25.00

X


9.00

11.24 ±0.25

63.00 ±0.25

85.48 ±0.25

X






25.00

25.00

Y

31.00 ±0.25

6.00 ±0.25



60.10

PLATES STACKED

All dimensions in mm

SECTION X-X