




NOTES:

1. MATERIAL: ANODIZED ALUMINUM PLATE WITH KAPTON LINER
2. UPPER LEFT CORNER NOTCHED FOR QUICK A1 IDENTIFICATION
3. APPROX. VOLUME: 400  $\mu$ L/WELL
4. PLATE WILL BE ENGRAVED WITH ALPHANUMERICS

REVISION HISTORY			
REV	DESCRIPTION	DATE	APPROVED
1A	INITIAL RELEASE	4/22/2026	ENGINEERING



DRAWN MD	4/22/2026	<b>ANALYTICAL</b> SALES AND SERVICES, INC © 2026	
email engineering@analytical-sales.com			
		TITLE	
		968040_rev1a_PUBLIC	
		SIZE B	REV 1A
		SCALE 1 : 1	SHEET 1 OF 1

P