



NOTES:
1. MATERIAL: GLASS, CLEAR
2. VOLUME: 2mL NOMINAL


REVISION HISTORY			
REV	DESCRIPTION	DATE	APPROVED
1A	INITIAL RELEASE	4/8/2022	ENGINEERING



A



A



SECTION A-A

P

DRAWN Matt Dabrowski email engineering@analytical-sales.com	4/8/2022	ANALYTICAL SALES AND SERVICES, INC © 2022		
		TITLE		
		96F620_2mL_flat_bot_insert_rev1a_P		
		SIZE B		REV 1A
		SCALE 3 : 1		SHEET 1 OF 1