


NOTES:

- 1. ELECTRICAL RATING:
VOLTS: 24
WATTS: 24
- 2. OPERATING TEMPERATURE RANGE:
5°C ABOVE AMBIENT TO 85°C

REVISION HISTORY			
REV	DESCRIPTION	DATE	APPROVED
1	INITIAL PUBLIC RELEASE	11/30/2023	ENGINEERING



END VIEW
SCALE 2 : 1



FRONT VIEW
SCALE 1 : 1

P

DRAWN Sebastian Sosa email engineering@analytical-sales.com	11/30/2023	ANALYTICAL SALES & SERVICES, INC © 2023	
DIMS IN INCH [MM]	TITLE HSI-10L_6in_Hot_Sleeve_rev1_PUBLIC		
	SIZE B		REV 1
	SCALE 3 : 4	SHEET 1 OF 1	