

Lumidox 96-Well Green LED Array




LUM96G

Optical Data per LED:

Dominant Wavelength..... 527 nm
Activation Current..... 1.5 mA
Maximum Radiant Flux..... 6.4 mW at 30 mA
Radiant Flux Reading..... 3.3 mW at 15 mA



Measured Radiant Flux vs. Current per LED



NOTE ¹: All information presented on this data sheet is given as measured per LED

NOTE ²: Measurements taken from the top of the domed cover

Maximum Radiant Flux in 96 Array: 616.3 mW