

# 2mL TrueTaper® 96-Well Collection Plate with 100uL Tapered-Bottom Wells

PN: 968820

## Dimensions (mm)

### Plate Size

X: 127.8 mm

Y: 85.5 mm

Height: 44 mm

### Top Left Well Location

X: 14.38 mm

Y: 11.24 mm

### Well Size

Square: 8.20 mm

Depth: 39.9 mm

### Rows

Number: 8

Spacing: 9 mm

### Columns

Number: 12

Spacing: 9 mm