

2.0mL 96-Well Flexi-Tier™ Block System

Opaque and Clear

PN: 96FP20 and 96FP20-C

Dimensions (mm)

Plate Size

X: 126.9 mm

Y: 84.9 mm

Height: 61.1 mm

Top Left Well Location

X: 13.45 mm

Y: 10.71 mm

Well Size

Diameter: 8.97 mm

Depth: 54.5 mm

Rows

Number: 8

Spacing: 9 mm

Columns

Number: 12

Spacing: 9 mm